

Next Generation Sensors for Pharmaceutical Process Analysis

*Ray D. Reid, R. Bhartia, Q. Nguyen, M. Reid, K. Sijapati, K. Nguyen, &
William Hug*

Photon Systems, Inc.

*FACSS/SciX
Providence RI
September 26, 2021*



1512 Industrial Park St., Covina, CA 91722 T: 626 967-6431 F: 626 967-5813

www.photonsystems.com

Outline

- ❑ Detection goals & methods for RCV and real-time continuous manufacturing
- ❑ Advantages of deep UV Raman & fluorescence detection
- ❑ Detection examples for pharma products
- ❑ Deep UV Raman & fluorescence instruments
- ❑ Chemical printer for NIST traceable chemical concentration calibration

Advantages of deep UV Raman & fluorescence detection



PHOTON
S Y S T E M S

1512 Industrial Park St., Covina, CA 91722 T: 626 967-6431 F: 626 967-5813

www.photonsystems.com



Advantages of Deep UV Detection vs Visible or IR?

- ❑ Non-contact, reagentless, no sample handling or preparation
- ❑ Excitation below 250 nm separates Raman & fluorescence spectral regions to enable
 - ✓ Clear Raman spectra with no obscuration or alteration by native fluorescence
 - ✓ No alteration of the fluorescence spectra by major Raman bands
 - ✓ The ability to simultaneously detect Raman and native fluorescence
- ❑ Much higher Raman sensitivity due to Rayleigh law and resonance Raman enhancement effects
- ❑ Fluorescence detection alone has much higher specificity when excitation is below 250 nm
- ❑ Detection of concentration of pharma materials in the low ng/cm² has been demonstrated
- ❑ Detection is solar blind, enabling detection in full daylight without interferences



Deep UV Raman Advantages


Primary advantages



- Spectral separation of Raman and fluorescence
- Increased sensitivity to UV-absorbing compounds (primarily aromatics, but others exhibit pre-resonance)
- Simplification of Raman spectra

Why Deep UV below 250nm?


When excitation < 250nm Raman and fluorescence spectral regions are separated



Sensitivity to Excitation Wavelength

Raman Spectra with Excitation at 248 nm versus 262 nm


(Example is G Agents)



Adapted from Christesen, SD et al. Appl Spec. 2008 Oct; 62(10):1078-83

Deep UV Fluorescence Spectra of 52 Compounds


with no baseline subtraction or compensation, Ex=248 nm



Chemical Differentiability

Using Deep UV Excited Fluorescence Alone

A single deep UV laser pulse determines the location of an unknown substance in this chemometric space



Combining the Sensitivity of Fluorescence & specificity of Raman

- ❑ Fluorescence is the most sensitive method of detection, over 10^6 to 10^8 times more sensitive than Raman, providing longer standoff distances and/or detection at lower concentrations
- ❑ Raman provides information about chemical bonds and functional groups, including those that do not fluoresce (aliphatics and simple compounds)
- ❑ Fluorescence data provides information about the overall electronic structure of target & substrate components (aromatics, ketones, aldehydes)

Raman Active		Weak Fluorescence		Strong Fluorescence	
Water	HMX	TDG	DMMP	C4	Microbes
Amino Acids	PETN	DIMP	TEPO	Semtex	Toxins/Proteins
Alcohols	RDX	Ammonia Nitrate		ANFOs	Narcotics
Aliphatics		Urea Nitrate			
DNA/RNA	TNT	Nitroglycerin		Aromatic Amino Acids	
Lipids	Perchlorates	Ketones/Aldehydes		PETN	VOCs


Detection Examples for Pharm Products

Pharma Applications and Products


TraC: Small Deep UV multi-channel fluorescence with high sensitivity



- Rapid Cleaning Verification RCV tool
- Feed Frame detector: low API concentrations high potency drug mixtures (in development)





RPL 200: Portable Deep UV high resolution Raman and fluorescence spectral analysis




- Enables Raman without fluorescence obscuration.
- Enables detection not possible with 785 or 1064 systems
- Microbial detection

Benlyn: dextromethorphan hydrobromide




248.6 nm data: Raw, unprocessed data: No Fluorescence

ChemCal: NIST Traceable chemical printer




- Enables printing of a priori concentrations of Pharma materials.
- Enables direct concentration calibration of many RCV tools
- Training and evaluation of Swabbing techniques.






PHOTON
S Y S T E M S

Raman/PL 200 Flow Cell Module




Modular front-end design to
fit a variety of applications

Flow cell module currently
being used for in-line waste
water analysis (NO₃,
Organics, etc.)



Surface analysis
module being used for
microbial detection



DUV Raman/PL 200


Features of the Deep UV Raman PL 200

- ❑ A deep UV Raman and fluorescence spectrometer
- ❑ With either computer-controlled stage for mapping or liquid flow cell for continuous manufacturing quality control
- ❑ Intended for OEM applications with dramatically smaller SWAP/C than other deep UV instruments on the market
- ❑ Avoids fluorescence interference or obscuration of Raman spectra
- ❑ Enables detection and quantification of Raman bands for a wide range of pharma ingredients not possible with 785 nm or 1064 nm Raman systems due to fluorescence
- ❑ SWAP: 18 x 20 x 42 cm, 10 kg, 60 W max (100-260VAC)



Raman PL 200

with various types of flow cells or cuvette holders





PHOTON
S Y S T E M S


1512 Industrial Park St., Covina, CA 91722 T: 626 967-6431 F: 626 967-5813

www.photonsystems.com

OTC Benylin: dextromethorphan hydrobromide




Benylin (white)
Sodium benzoate (green)
Glycerol (red)




OTC Children's Motrin (ibuprofen)–Bubblegum Flavor

Ex = 248 nm

Children's Motrin-Bubblegum flavor




Pure Ibuprofen



OTC Children's Tylenol (acetaminophen) w Various Flavors

Ex = 248.6 nm Raw results. No baseline compensation.



Deep UV Trace Chemical (TraC) Sensor

LODs < 1 $\mu\text{g}/\text{cm}^2$ Wt. 1.5 lbs



PHOTON
S Y S T E M S

1512 Industrial Park St., Covina, CA 91722 T: 626 967-6431 F: 626 967-5813 www.photonsystems.com

Features of Deep UV TraC RCV Sensor

- ❑ Fully self-contained RCV sensor with embedded microprocessor for instrument control, data processing, real-time data storage, and display
- ❑ Able to measure trace concentrations on curved surfaces, corners, crevices, screens, grates
- ❑ Sample rate > 10 samples/s with time-stamped real-time recording
- ❑ Hi sensitivity: able to detect concentrations of APIs < 1 µg/cm²
- ❑ Large working distance: 0 to 2 cm
- ❑ Sampling area: 0.25 cm²
- ❑ Non-contact sensing with large working distance (0 to 2 cm)
- ❑ Hand-Held: < 0.7 kg (1.5 lbs)
- ❑ Small: 7.6 x 8.9 x 19 cm
- ❑ Long battery lifetime: > 40 hours full power; > 120 hours standby
- ❑ Startup time < 10 s
- ❑ GMP & Intrinsically safe



PHOTON
S Y S T E M S

1512 Industrial Park St., Covina, CA 91722 T: 626 967-6431 F: 626 967-5813

www.photonsystems.com

Quality Control in Continuous Manufacturing

The goal

Provide instrumentation for real-time detection of the key ingredients during continuous flow manufacturing

The solution

A miniature deep UV instrument that enables quantification of API extending the current NRI levels of ~2% volume to well below .1% mixture. Feed Frame



PHOTON
S Y S T E M S

1512 Industrial Park St., Covina, CA 91722 T: 626 967-6431 F: 626 967-5813

www.photonsystems.com

Detector Choices for Moving Powder

Deep UV application Feed Frame

- **Problem:**

- NIR and FITR techniques are limited in most powder mixtures to > 2% bulk ratios. Absorption techniques and excipient interferences limits dynamic range.

- **Goal:**


- Explore other techniques focused in Deep UV spectroscopy (fluorescence and Raman) to extend in process control of high potency drugs to better than 0.1% bulk ratios.

- **Solution:**

- A handheld size, deep UV Raman/Fluorescence instrument, that avoids spectral obscuration enabling the advantage of both spectroscopic techniques.



Fluorescence

- High sensitivity with API
- Typically, 6 order magnitude more sensitivity than Raman
- Easy to configure fluorescence detection to optimize sensitivity and specificity
- Meets GMP requirements



Feed Frame
instrument under
development

API #1 Fluorescence




API #1

LOD calculation is 3*stdev signal to background 0%API

LOQ is 10*stdev signal to background 0%API

Note stdev for channel 4 is quite small resulting in very low LOD/LOQ

API #1	(Placebo)	TraC300		
Pulse #	CH1	CH2	CH3	CH4
1	1158016	4330302	15426030	562882
2	1158790	4324586	15414404	563500
3	1154584	4324712	15407616	566510
4	1158775	4313648	15405742	563924
5	1150324	4314752	15389699	555454
6	1162153	4310587	15384947	557960
7	1150843	4305384	15379319	558977
8	1153033	4306050	15377898	562912
9	1153379	4299269	15359867	567778
10	1151667	4301042	15346608	559587
avg	1155156	4313033	15389213	561948
stdev	4022	10635	24772	3879
stdev*3	12065	31905	74315	11636
stdev*10	40216	106349	247716	38786
graph equation				graph equation
x^2				-12049151022
x				650129212.8
c				600259.64
LOD in % using -				5.393853%
LOQ in % using -				5.389671%
LOD in % using +				0.001790%
LOQ in % using +				0.005973%







API #1

TraC300 Repeatability on API #1 samples

API #1 TraC-300 Repeatability w/ 15% Error Bars



API #1 Repeatability %from AVG




- TraC set stationary on the ChemCal and the samples were manually moved underneath the TraC nose.
- Repeatability of the various concentrations range from 3.0% to 6.7%, which are well within the $\pm 15\%$ repeatability limits previously set for API TraC development.

TraC: Saccharin concentration curves

Fluorescence response Feed Frame application

TraC

- 3.00% saccharin in 5 standard excipients is 100% mixture
- Dilution from 125% (3.75w/w) to 3% (0.1%w/w)
- We would estimate the LOD with SN of 3 to be ~0.04% W/w Saccharin.
- Roughly 20 times more sensitive than NIR techniques currently used.
- Currently under study with ETC (Emerging Technology Consortium) a Global Pharma group



Rapid Cleaning Verification

The goal

To augment or replace the present swab & test method for equipment cleaning verification with a faster and better controlled and documented method.

The solution


A handheld device that quantifies trace amounts of API in real time on manufacturing surfaces.

Result: Significantly reducing production down time.




Typical Concentration Calibration Curve & Time Stability

Ibuprofen Chemical Concentration Calibration



Ibuprofen Signal Time Study



Why do you need a chemical printer?

- ❑ Test/calibrate/validate future RCV tools using NIST traceable method
- ❑ Test/calibrate/validate CURRENT cleaning tools & methods
 - ✓ Create concentration curves for swabbing.
 - ✓ Test/train swabbing personnel with accurate areal concentrations.
 - ✓ Test recovery from various surfaces/topologies with different swabbing media.
- ❑ Create coupons for visual /hotspot detection of API.
 - ✓ Hotspot detection.
 - ✓ Train personnel on visual inspection limits.
- ❑ Perform all of these with single or multiple chemicals on a single coupon or coupons.
 - ✓ Detergent + API.
 - ✓ Excipient + API.
 - ✓ Detergent + API + Excipient.




PHOTON
S Y S T E M S

1512 Industrial Park St., Covina, CA 91722 T: 626 967-6431 F: 626 967-5813

www.photonsystems.com


Chemical Calibration Printing



NIST Traceable Chemical Concentration Calibration

ChemCal: A chemical printer, mapper, & calibrator

- Creates up to 16 coupons with *a priori* known, NIST traceable, concentrations of many different chemicals, including APIs, detergents, excipients, etc. on Pharma-type surfaces or quartz crystal microbalance elements for the purpose of performing calibration of hand-held trace chemical sensors for rapid cleaning validation.
- Prints and detects on curved or flat surfaces, corners, grates, screens, etc.



Coupons

TraC sensor

ChemPrint head




Operational Scenario:

- ✓ Load APIs, etc in Eppendorf rack. Up to 21.
- ✓ Load coupons onto tray.
- ✓ Press Start.
- ✓ The system prints, scans, & outputs a full calibration curve in under 3 hours. (Prints and scans 16 coupons.)

Printing


Up to 5x5 cm grid

NIST Tracible Syringe



High surface tension solvent. Individual drops.

TraC Beam Diameter



Low surface tension solvent "Fish Scale Pattern" Overlapping drops



PHOTON
S Y S T E M S

1512 Industrial Park St., Covina, CA 91722 T: 626 967-6431 F: 626 967-5813

www.photonsystems.com

Summary

- ❑ Several opportunities exist for pharma and other chemical manufacturing RCV and continuous flow manufacturing instrumentation using deep UV excitation below 250 nm.
- ❑ Excitation below 250 nm provides fluorescence-free Raman and Raman-free fluorescence detection simultaneously, enabling both modes of detection to provide more accurate information about a trace substance on a surface or in a liquid.
- ❑ Combined Raman & fluorescence detection method enhance both sensitivity and specificity in identifying unknown targets
- ❑ Detection of Raman & fluorescence in the deep UV can be accomplished using low energy lasers without major alteration or damage/ignition of targets.
- ❑ Surface detectors need a method of accurate chemical concentration calibration, common to all methods of surface detection. We call this instrument **ChemCal**.



Questions ?



1512 Industrial Park St., Covina, CA 91722 T: 626 967-6431 F: 626 967-5813 www.photonsystems.com